

New Reservoir developed in pursuant of utmost usability.


PP PIPETTING RESERVOIR

20mL

Maximum volume
up to 50mL

NEW "Liquid Sweeper" has realized comfortable pipetting, developed focusing on usability and a feature includes the thought-through angle for gathering a liquid efficiently.

Liquid Sweeper




8 grooves, a very unique feature called "Liquid Sweeper", transfer a liquid drop from the top of the tip quickly.



Drop-on Cover



The most suitable angle for tip touch



The calculated slope angle gathers a liquid almost without its dead volume.

Convenient graduation



Easy-to-wash round corners



Autoclavable

Note: Minor deformation may occur due to thermal contraction, but it does not affect usage.

Wide size to cover 12ch pipette

Cat. No.	Item	Unit	Case
1790-050-30	20mL, PP pipetting reservoir	30pcs.	10units
1791-050-S	20mL, PP pipetting reservoir, individual package (with cover)	Sterilized 10pcs.	5units

Something Different.

WATSON BIO LAB
MADE IN JAPAN SINCE 1988

<https://watsonbiolab.com>

Manufactured and sold by:

FUKAEKASEI CO., LTD.

Head Office: 2-2-7 Murotani, Nishi-ku, Kobe 651-2241, JAPAN TEL +81-78-991-4477 FAX +81-78-991-4491
E-mail: info@watson.co.jp